













































Ordering Information

Curve	Tip	5 Fr	6 Fr	7 Fr	8 Fr
 Extra Back-up Left	XB3.0	53510004	53610001	53710004	53810004
 Extra Back-up Left	XB3.5	53510008	53610002	53710008	53810008
 Extra Back-up Left	XB4.0	53510012	53610003	53710012	53810012
 Judkins Left	JL3.0	53510016	53610022	53710016	53810016
 Judkins Left	JL3.5	53510020	53610004	53710020	53810020
 Judkins Left	JL4.0	53510024	53610005	53710024	53810024
 Judkins Left	JL4.5	53510028	53610006	53710028	53810028
 Judkins Left	JL5.0	53510032	53610007	53710032	53810032
 Judkins Right	JR3.0	53510036	53610023	53710036	53810036
 Judkins Right	JR3.5	53510040	53610008	53710040	53810040
 Judkins Right	JR4.0	53510044	53610009	53710044	53810044
 Judkins Right	JR4.5	53510048	53610010	53710048	53810048
 Judkins Right	JR5.0	53510052	53610011	53710052	53810052
 Amplatz Left	AL1	53510056	53610012	53710056	53810056
 Amplatz Left	AL2	53510060	53610013	53710060	53810060
 Amplatz Left	AL3	53510064	53610014	53710064	53810064
 Amplatz Right	AR	53510068	53610015	53710068	53810068
 Amplatz Right	AR1	53510072	53610016	53710072	53810072
 Amplatz Right	AR2	53510076	53610017	53710076	53810076
 Williams	3DRC	53510080	53610018	53710080	53810080
 Multipurpose	MPA	53510084	53610019	53710084	53810084
 Multipurpose	MPB	53510088	53610021	53710088	53810088

Curve (with side hole)	Tip	5 Fr	6 Fr	7 Fr	8 Fr
 Extra Back-up Left	XB3.0	53510002	53610026	53710002	53810002
 Extra Back-up Left	XB3.5	53510006	53610028	53710006	53810006
 Extra Back-up Left	XB4.0	53510010	53610031	53710010	53810010
 Judkins Left	JL3.0	53510014	53610034	53710014	53810014
 Judkins Left	JL3.5	53510018	53610037	53710018	53810018
 Judkins Left	JL4.0	53510022	53610040	53710022	53810022
 Judkins Left	JL4.5	53510026	53610043	53710026	53810026
 Judkins Left	JL5.0	53510030	53610046	53710030	53810030
 Judkins Right	JR3.0	53510034	53610049	53710034	53810034
 Judkins Right	JR3.5	53510038	53610052	53710038	53810038
 Judkins Right	JR4.0	53510042	53610055	53710042	53810042
 Judkins Right	JR4.5	53510046	53610058	53710046	53810046
 Judkins Right	JR5.0	53510050	53610061	53710050	53810050
 Amplatz Left	AL1	53510054	53610064	53710054	53810054
 Amplatz Left	AL2	53510058	53610067	53710058	53810058
 Amplatz Left	AL3	53510062	53610070	53710062	53810062
 Amplatz Right	AR	53510066	53610073	53710066	53810066
 Amplatz Right	AR1	53510070	53610076	53710070	53810070
 Amplatz Right	AR2	53510074	53610079	53710074	53810074
 Williams	3DRC	53510078	53610082	53710078	53810078
 Multipurpose	MPA	53510082	53610085	53710082	53810082
 Multipurpose	MPB	53510086	53610088	53710086	53810086

CE 1639
APT Medical Inc.

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. of China
 Post code: 518057 Email: international@aptmed.com
 Tel: +86 755-83591013 Fax: +86 755-83480508 en.aptmed.com
 PN: ENG-March Guiding Catheter-210X285X4P-20201030

March[®] Guiding Catheter



➤ Excellent Kink Resistance

Polymer material with unique braid design, provides excellent kink resistance.




March®

Brand M


➤ Superior Support

Multi-density braid design provides superior back-up support.




➤ Radiopaque Soft Tip

Radiopaque soft tip ensures atraumatic navigation.



➤ Optimal Flexibility

Torque Response Comparison




1:1 torque response promotes quicker catheter seating with less whipping.

Source: internal bench test data comparisons.

➤ Large Lumen Design

Provides smooth and compatible pathway for multiple and various devices.



	5 Fr	6 Fr	7 Fr	8 Fr
I.D. (Inch/mm)	0.058/1.47	0.071/1.80	0.081/2.06	0.093/2.36
O.D. (Inch/mm)	0.069/1.75	0.0815/2.07	0.093/2.36	0.105/2.67