























Ordering Information

Guiding Catheter

Curve	Tip Code	5.5 F	6.5 F	5.5 F (Side Hole)	6.5 F (Side Hole)
 Extra Back-up Left	XB3.0	53110004	53210004	53110002	53210002
 Extra Back-up Left	XB3.5	53110008	53210008	53110006	53210006
 Extra Back-up Left	XB4.0	53110012	53210012	53110010	53210010
 Judkins Left	JL3.0	53110016	53210016	53110014	53210014
 Judkins Left	JL3.5	53110020	53210020	53110018	53210018
 Judkins Left	JL4.0	53110024	53210024	53110022	53210022
 Judkins Left	JL4.5	53110028	53210028	53110026	53210026
 Judkins Left	JL5.0	53110032	53210032	53110030	53210030
 Judkins Right	JR3.0	53110036	53210036	53110034	53210034
 Judkins Right	JR3.5	53110040	53210040	53110038	53210038
 Judkins Right	JR4.0	53110044	53210044	53110042	53210042
 Judkins Right	JR4.5	53110048	53210048	53110046	53210046
 Judkins Right	JR5.0	53110052	53210052	53110050	53210050
 Amplatz Left	AL1	53110056	53210056	53110054	53210054
 Amplatz Left	AL2	53110060	53210060	53110058	53210058
 Amplatz Left	AL3	53110064	53210064	53110062	53210062
 Amplatz Right	AR	53110068	53210068	53110066	53210066
 Amplatz Right	AR1	53110072	53210072	53110070	53210070
 Amplatz Right	AR2	53110076	53210076	53110074	53210074
 Williams	3DRC	53110080	53210080	53110078	53210078
 Multipurpose	MPA	53110084	53210084	53110082	53210082
 Multipurpose	MPB	53110088	53210088	53110086	53210086

Hemostasis Introducer

French Size	Introducer needle	Scalpel	Sheath Introducer				Marker Tip Sheath Introducer	
			7cm	11cm	15cm	Guidewire	25cm	Guidewire
5.5F	22G 45mm	11 P	98050071	98050111	98050151	Straight 0.025" 45cm	98050251	Straight 0.025" 80cm
6.5F	22G 45mm	11 P	98060071	98060111	98060151	Straight 0.025" 45cm	98060251	Straight 0.025" 80cm


French Size	Introducer needle	Guidewire		Sheath Introducer without Scalpel			Sheath Introducer with Scalpel		
		Tip	Diameter	Length	11cm	15cm	Scalpel	15cm	Scalpel
5.5F	21G 70mm	Straight	0.018"	45cm	97050111	97050153	—	97050151	11 P
6.5F	21G 70mm	Straight	0.018"	45cm	97060111	97060153	—	97060151	11 P

French Size	Introducer needle	Guidewire		Sheath Introducer without Scalpel			Sheath Introducer with Scalpel		
		Tip	Diameter	Length	11cm	15cm	Scalpel	15cm	Scalpel
5.5F	21G 40mm	Straight	0.018"	45cm	97050113	97050157	—	97050155	11 P
6.5F	21G 40mm	Straight	0.018"	45cm	97060113	97060157	—	97060155	11 P

CE 1639



6.5F Hemostasis introducer and guiding catheter




APT Medical Inc.

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. of China
 Post code: 518057 Email: international@aptmed.com
 Tel: +86 755-83591013 Fax: +86 755-83480508 en.aptmed.com
 PN: ENG-6.5F Hemostasis introducer and guiding catheter-210X285X4P-20201030

➤ 6.5F Hemostasis Introducer Thin Wall Design

- Thin wall design, equivalent outer diameter between 6.5F APT hemostasis introducer and conventional 6F introducer.

Comparison of hemostasis introducer wall



➤ 6.5F Guiding Catheter 0.075" Large Lumen

- Optimised solution for transradial intervention.
- Compatible with more interventional devices, including atherectomy with bigger rotablater, kissing technique with double balloon or stents, meet variety of clinical needs and challenge from complex lesions.


Compatible	6.5F APT hemostasis introducer and guiding catheter	Conventional 6F introducer and guiding catheter
Balloon & Stent	+	+
Stent & Stent	+	±
IVUS Catheter	+	+
Thromboectomy Catheter	+	+
Rotablater	≤1.75mm	<1.5mm



6.5F APT introducer and guiding catheter



6F conventional introducer and guiding catheter



Half French larger with big difference for you and your patients