


# Ordering Information

## Transseptal Sheath


Frenchsize	Curve	Sheath Length (cm)	Dilator Length (cm)	Guidewire (0.025" 150cm 3mm J curve)		
				Spring	Hydrophilic	PTFE
6F	T0	45	52	/	94080450	94080454
6F	T1	45	52	/	94060451	94060455
6F	T2	45	52	94060453	94060452	94060456
6F	T3	45	52	/	94060453	94060457
7F	T2	45	52	94070453	94070452	94070454
7F	T3	45	52	/	94070451	94070455
8F	T2	45	52	94080453	94080452	94080454
8F	T3	45	52	/	94080451	94080455

Frenchsize	Curve	Sheath Length (cm)	Dilator Length (cm)	Guidewire (0.032" 150cm/180cm 3mm J curve)		
				Spring	Hydrophilic	PTFE
6F	T2	60	67	94060603	94060602	94060604
6F	T3	60	67	/	94060601	94060605
7F	T2	60	67	94070603	94070602	94070604
7F	T3	60	67	/	94070601	94070605
8F	T0	62	67	94080622	94080620	94080627
8F	T0	80	85	94080802	94080800	94080809
8F	T1	62	67	94080621	94080623	94080628
8F	T1	80	85	94080801	94080803	94080808
8F	T2	60	67	94080603	94080602	94080604
8F	T3	60	67	/	94080601	94080605
8F	T3	62	67	/	94080624	94080629
8F	T3	80	85	/	94080804	94080807
8.5F	T0	62	67	94080626	94085620	94085624
8.5F	T0	80	85	94080806	94085800	94085804
8.5F	T1	62	67	94080625	94085621	94085625
8.5F	T1	80	85	94080805	94085801	94085805
8.5F	T3	62	67	/	94085623	94085626
8.5F	T3	80	85	/	94085803	94085806
10F	T1	60	67	/	94100601	94100605
10F	T2	60	67	94100603	94100602	94100607
10F	T3	60	67	/	94100604	94100606
12F	T1	60	67	/	94120601	94120605
12F	T2	60	67	94120603	94120602	94120607
12F	T3	60	67	/	94120604	94120606

Guidewire length: 150cm    Guidewire length: 180cm



3D curve



## Transseptal Needle

REF	Model	Usable Length(mm)	Stylet Length(mm)	Curve
94180710	18G	710	765	
94180890	18G	890	920	
94190560	19G	560	590	
94180711	18G	710	765	
94180891	18G	890	920	
94190561	19G	560	590	

Please contact us for more information:




### APT Medical Inc.

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. of China

Post code: 518057 Email: international@aptmed.com


Tel: +86 755-83591013 Fax: +86 755-83480508 en.aptmed.com

PN: ENG-Transseptal needle and sheath-210X285X4P-20201030



# Transseptal Needle and Sheath

# Transseptal Needle and Sheath



## Product Benefits

- ⊙ Nylon elastic material of sheath ensures good curve within the heart
- ⊙ PTFE inner layer of sheath enhances smooth passage of catheter
- ⊙ Special plastic injection moulding handle maintains it's reliability

